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# Samsung Health Monitor

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The Product described in below is medical device

# Samsung Health Monitor

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### (ECG Monitor App)

129, Samsung-ro,  
Yeongtong-gu,  
Suwon-si, Gyeonggi-do,  
Republic of Korea

**Samsung Electronic Co., Ltd.**

[www.samsung.com](http://www.samsung.com)

Please read this instruction  
for use carefully to use the  
product properly

## Instruction for Use

2020. 06. 25 (Rev. A)

## 1 Intended use

The Samsung Health Monitor ECG Monitor(ECG app) labels atrial fibrillation (AFib) or sinus rhythm by using ECG sensor on mobile platform (Smart Watch) and the app is intended to create, record, store, transfer, and display a single-channel electrocardiogram (ECG), similar to a Lead I ECG, for informational use only in adults 22 years and older.

## 2 Product overview

The Samsung ECG app is an over-the-counter (OTC) software-only, mobile medical application operating on a compatible Samsung Galaxy Watch and Phone. The app is intended to create, record, store, transfer, and display a single-channel electrocardiogram (ECG), similar to a Lead I ECG, for informational use only in adults 22 years and older. Classifiable traces are labeled by the app as either atrial fibrillation (AFib) or sinus rhythm with the intention of aiding heart rhythm identification; it is not intended to replace traditional methods of diagnosis or treatment. The app is not intended for users with other known arrhythmias and users should not interpret or take clinical action based on the device output without consultation of a qualified healthcare professional.

## 3 Warnings and Cautions

### WARNINGS

Please follow the warnings listed below.

	The ECG Monitor App cannot diagnose cardiac conditions or look for signs of a heart attack.
	The ECG Monitor App is not meant to replace traditional methods of diagnosis or treatment by a qualified healthcare professional. If you think you are having a medical emergency, contact your local emergency services immediately.
	DO NOT change your medications or dosage without first consulting your doctor.
	ECG rhythm-classification results displayed by the ECG Monitor App are for <b>informational purposes only</b> . DO NOT interpret or take clinical action based on the ECG App measurements without first consulting with a qualified healthcare professional.
	DO NOT use the ECG Monitor App if you have an implanted pacemaker, implanted cardiac defibrillator, or other implanted electronic devices.
	DO NOT use the ECG Monitor App if you have known arrhythmia other than Atrial Fibrillation.

**CAUTIONS**

Please follow the precautions below. If you do not follow them, the ECG App may not be able to record accurate measurements.

	DO NOT take ECG measurements during any physical activity.
	DO NOT use the ECG App if you are younger than 22 years old.
	DO NOT take ECG measurements when the Galaxy Watch is close to strong electromagnetic fields (for example, magnetic resonance imaging (MRI) or X-Ray equipment, electromagnetic anti-theft systems, and metal detectors).
	For security reasons, always pair the Galaxy phone and Galaxy watch via Bluetooth in a private (home-based) setting. It is NOT recommended that pairing be done in a public space
	DO NOT take ECG measurements during a medical procedure (for example, surgery or external defibrillation procedures).
	DO NOT take ECG measurements when the Galaxy Watch is outside of the following conditions: <ul style="list-style-type: none"> <li>• Temperature: 54 °F – 104 °F/12 °C – 40 °C</li> <li>• Humidity: 30% – 90% relative humidity</li> </ul>
	Make sure that the finger used to record the ECG is uncovered and free of any bruises, scars, or cuts.
	Certain physiological conditions can prevent some people from having a strong enough signal for the ECG App to detect and analyze.

**4 Getting Started**

**Note:** The ECG Monitor App requires a Samsung Galaxy Watch paired with a compatible Samsung Galaxy phone (check <https://www.samsung.com/sec/apps/samsung-health-monitor> for a list of compatible devices).

All Samsung Galaxy mobile devices with Android N and above and Samsung Galaxy Active 2 with Firmware 4.0.0.8 and later or later Samsung Galaxy Watch models are supported

1. Make sure that your Galaxy Watch is paired with your Galaxy phone using Bluetooth.
2. If the ECG Monitor App is installed on your Galaxy phone, open the Samsung Health Monitor app and follow the onscreen instructions to set up your profile.

If the ECG Monitor App is not installed on your Galaxy phone, install the Samsung Health Monitor App from the Galaxy store and then install the latest update for the Galaxy Wearable software. After installing the update, open the Samsung Health Monitor App and follow the onscreen set-up instructions.

## 5 Recording an ECG Using Your Galaxy Watch

---

1. Open the ECG Monitor app on your Galaxy Watch.
2. Make sure the Galaxy Watch is snug on your wrist.
3. Rest your forearms comfortably on a table and then rest a fingertip of the opposite hand *lightly* on the Galaxy Watch's Home key for 30 seconds. Remain still and do not talk while the Galaxy Watch takes the measurement. When the recording process has finished, the ECG results will appear on your watch. Note: **Don't Press** the Home key during recording; pressing the Home key during recording will stop the recording.
4. Scroll up and down to see your ECG result. If you are experiencing symptoms, tap **Add** and on the **Symptoms** screen, scroll through the list and tap each symptom you are experiencing. Then tap **Save** to add the symptoms to your ECG report and to return to the **Results** screen.
5. Tap **Done** to go back to the **Start** screen. You can review your ECG report in the ECG Monitor app on your phone. (Note that syncing requires that the watch and phone are paired using Bluetooth.)

## 6 Reviewing Your ECG Reports on Your Galaxy Phone

---

After you take an ECG, the ECG data is synced to the paired Galaxy phone where the ECG App analyzes the data and creates a PDF report.

1. Open the Health Monitor app on your Galaxy phone.
2. If Health Monitor App menu is available, Select the ECG Monitor tab in the Health Monitor app menu.
3. Tap the thumbnail to see your latest ECG report, or tap **View history** to see all the available reports.

You may get any of the four following results:

- **Sinus Rhythm**—This result means that during the recording the heart was beating in a regular rhythm with a heart rate of 50–100 beats per minute (BPM).
  - **Caution:** A Sinus Rhythm result does not guarantee that you are not experiencing an arrhythmia or other health condition. If you are not feeling well, contact your doctor
- **Atrial Fibrillation (AFib)**—This result means that during the recording the heart was beating in an irregular rhythm with a heart rate of 50–120 BPM. If you get this result, contact your doctor for guidance.
  - **Note:** Because the occurrence of asymptomatic atrial fibrillation is not constant and can fluctuate daily, periodic reevaluation is necessary. (Source: 2021 Korean Heart Rhythm Society Guidelines for Screening and Management of Subclinical Atrial Fibrillation)

- **Inconclusive**—This result means that the ECG recording could not be classified because the heart rate was either too high or too low, or the rhythm was not Atrial Fibrillation or Sinus Rhythm. If you get this result repeatedly, contact your doctor.
  - **Conditions that lead to an Inconclusive result**
    - The heart rate during the recording was less than 50 BPM, *or*
    - the heart rate was greater than 100 BPM and rhythm was not Atrial Fibrillation, *or*
    - the heart rhythm was not Sinus or Atrial Fibrillation, *or*
    - the heart rhythm was Atrial Fibrillation and the heart rate was greater than 120 BPM
  
- **Poor Recording**—This result means the ECG App was not able to analyze the data. A poor recording usually happens because your body moved during recording or the Galaxy Watch did not have enough skin contact with your wrist or finger. If you get this result repeatedly, refer to “Recording an ECG Using Your Galaxy Watch,” above

## 7 Safety and Performance

---

The clinical validation for ECG Monitor App tested its accuracy in detecting Sinus Rhythm and Atrial Fibrillation in ECG recordings for 544 subjects. The ECG Monitor App rhythm classifications were compared with rhythm classifications performed by a board-certified cardiologists using a 12-lead ECG. The ECG Monitor App had a sensitivity of 98.1% in detecting AFib and specificity of 100% in classifying sinus rhythm for all classifiable recordings.

During clinical trials, 16.8% of ECG recordings were either inconclusive or poor recording. If all of these recordings are included, the sensitivity of the ECG monitor App is 87.1% and specificity to detect sinus rhythm is 82.5%. Real-world performance may have more inconclusive and poor recordings.

The ECG PDF report was compared against a standard Lead 1 ECG for key intervals (PR, RR), QRS duration and amplitude to ensure its accuracy. No adverse events were reported during this clinical trial.

## 8 Troubleshooting

---

If you have problems using your ECG App, check the table below for possible solutions.

Problem	Solution
I cannot find the ECG Monitor App on my Galaxy Watch.	Make sure that your Galaxy Watch is compatible with the ECG Monitor App by checking the list on <a href="https://www.samsung.com/sec/apps/samsung-health-monitor">https://www.samsung.com/sec/apps/samsung-health-monitor</a> . If your Galaxy Watch is compatible, install the Samsung

	<p>Health Monitor App on the phone and upgrade your Galaxy Wearable software to the latest Maintenance Release (MR).</p>
<p>I cannot record my ECG</p>	<p>Dry or cool skin, not enough skin contact, hairy wrist, pressing the Home key or user movement can cause the ECG Monitor app to shut down, to not start recording, or to stop recording. Before taking an ECG:</p> <ul style="list-style-type: none"> <li>• Put some moisturizing lotion on your hands and wrist.</li> <li>• If your hands and arms are cold, remove the Galaxy Watch and rub the wrist gently to warm up the skin. Put the Galaxy Watch back on and try again to record an ECG.</li> <li>• Try wearing the watch on the underside of the wrist to create better skin contact.</li> </ul> <p>During the recording process make sure</p> <ul style="list-style-type: none"> <li>• Your forearms are resting on a table.</li> <li>• Your fingertip fully covers the Home key.</li> <li>• You do not move your arms or talk.</li> <li>• You do not press the Home key.</li> <li>• Your finger is resting <i>lightly</i> on the Home key for the 30 seconds it takes to do the recording.</li> </ul> <p>If the ECG App continues to shut down as soon as you open it, it may be a software issue. Restart the Galaxy Watch and try again. If you still cannot get a successful recording, it could be because the signal is not strong enough to record your ECG.</p>

<b>Problem</b>	<b>Solution</b>
<p>I don't see my ECG result in the ECG Monitor App.</p>	<p>Make sure your Galaxy Watch is paired via Bluetooth with the Galaxy phone using the Galaxy Wearable application. If you still don't see your results on the ECG Monitor App home screen, tap <b>Sync</b> in the upper right corner.</p>
<p>I am repeatedly getting a Poor Recording result.</p>	<p>A Poor Recording result may be due to "noise" in the signal caused by dirt, lack of moisture, insufficient skin contact, or user movement during the ECG measurement.</p> <ul style="list-style-type: none"> <li>• Clean the back of the watch and the fingertip you're using for the measurement.</li> <li>• Make sure the Galaxy Watch is snug on your wrist, rest your forearms comfortably on a table, and then take the ECG.</li> </ul>

	<ul style="list-style-type: none"> <li>You may need to wear the watch on your wrist slightly away from your wrist bone toward elbow to ensure good contact between your skin and watch.</li> <li>Certain physiological conditions can prevent some people from having a strong enough signal for the ECG Monitor App to detect and analyze.</li> </ul>
<p>I am repeatedly getting an Inconclusive result.</p>	<p>If you think you are having a medical emergency, contact your emergency services immediately.</p> <ul style="list-style-type: none"> <li>An Inconclusive result may be from a potentially high or low heart rate or another heart rhythm that is not AFib or Sinus Rhythm.             <ul style="list-style-type: none"> <li>Stay still for 5–10 minutes and try taking the ECG again if you did exercise, had food/beverage that includes caffeine recently or was in a stressful condition.</li> </ul> </li> <li>An Inconclusive result could be because the heart rate is more than 120 BPM and the rhythm is Atrial Fibrillation.</li> </ul> <p>If you get Inconclusive result repeatedly, contact your doctor.</p>
<p><b>Problem</b></p>	<p><b>Solution</b></p>
<p>My ECG waveform appears upside down.</p>	<p>An inverted waveform could be due to the wrong wrist selection. Please check your wrist setting and select the wrist that the ECG is recorded from.</p>
<p>Wrist orientation is reset to left</p>	<p>The wrist orientation is reset to the default setting (left) with Wear OS update. Please check your wrist setting and select the wrist that the ECG is recorded from.</p>